**Parvovirus B19
Immunofluorescent Assay**

An immunofluorescence assay for the detection of anti-Parvovirus B19 IgG and IgM antibodies.

Table of Contents

Intended Use	3	Specimen and Reagent Preparation	9
Introduction	3	Assay Procedure	10
Assay Principle	5	Interpretation of Results	11
Kit Components	5	Limitations of Use	13
Materials Provided	5	Performance Characteristics	14
Additional Materials Required	6	Summary of Parvovirus IFA Procedure	22
Precautions	7	Interpretation of symbols	22
Safety	7	Bibliography	165
Procedural	8		
Storage and Stability	9		
Specimen Collection and Storage	9		

Intended Use

The Parvovirus B19 Immunofluorescent Assay is intended for the qualitative and semi-quantitative detection of anti-Parvovirus B19 IgG and IgM antibodies in human serum. The test is indicated for use in all women where there is a suspicion of exposure to Parvovirus B19 to determine the presence of a recent or current infection.

Introduction

Parvovirus B19 was first identified as a human pathogen in 1975 and has subsequently been shown to be the causative agent of a number of clinical conditions such as rash, arthralgia and foetal damage.^{1,2,3} Parvovirus B19 infection in adults, especially women, may cause acute arthritis which can persist for some time.⁴ Infection can lead to life threatening anaemia in immunocompromised patients and individuals with underlying haemolytic disorders such as sickle cell disease.^{5,6} The virus is an icosahedral, non-enveloped virus of 18 – 25nm

diameter and comprises a linear single stranded DNA genome (5.5kb) which is encapsulated within an outer capsid.^{7,8} The viral capsid is composed of two structural proteins, namely VP1 (83kDA) and VP2 (53kDA). Parvovirus B19 infection is normally acquired by direct contact with respiratory secretions and normally occurs in localised outbreaks during the winter and spring months.⁸

It is now accepted that seronegative women are susceptible to Parvovirus B19 infection.^{9,10} The majority of pregnancies during which Parvovirus B19 infection occurs result in delivery of a healthy foetus at term.^{10,11,12} However, infection during pregnancy presents the risk of transmission to the foetus which may result in hydrops foetalis or intrauterine death. Estimates in the literature, for the rate of foetal death following maternal infection range between 1 and 9%.^{10,13,14} It has been suggested that because Parvovirus B19 replicates predominantly in red blood cell precursors, infection during pregnancy can lead to foetal death due to severe foetal anaemia. It is thought that this

severe anaemia, whereby haemoglobin levels fall to less than 2g/dl, is the primary cause of foetal hydrops.^{15,16}

The symptoms associated with Parvovirus B19 infection only become apparent after the viraemic (contagious) stage has terminated¹⁰. Furthermore, it is known that there is an increased risk of transmission in situations where close contact between individuals is likely, such as schools, day care centres and hospitals. Centres for Disease Control (CDC) do not recommend that persons exhibiting signs of Parvoviral infection (e.g. erythema infectiosum) be excluded from such environments. It is, however, recommended that all relevant individuals are made aware of the possibility of disease transmission.¹⁰

Consequently, it is important to identify the Parvovirus B19 antibody status in individuals who may be at risk of infection from, or who have been infected with, Parvovirus B19.

Serological assays are the mainstay of B19 diagnosis and Biotrin International offers a unique portfolio of Parvovirus B19 assays suitable for routine laboratory diagnosis.

- Parvovirus B19 IgG and IgM IFA
Catalogue No. V119IF
- Parvovirus B19 IgM EIA
Catalogue No. V619IM
- Parvovirus B19 IgM EIA (Automate)
Catalogue No. V619IMAUT
- Parvovirus B19 IgG EIA
Catalogue No. V519IG
- Parvovirus B19 IgG EIA (Automate)
Catalogue No. V519IGAUT

Assay Principle

The Biotrin Parvovirus B19 Immunofluorescent Assay utilises the indirect immunofluorescence antibody technique first described by Coons et al.¹⁷ Patient serum is incubated with Parvovirus B19 recombinant VP1 antigen in insect cells stabilised on a glass slide.^{18,19} If anti-Parvovirus B19 antibodies are present in the sample, a stable complex is formed with the antigen. Bound antibody is then reacted with a fluorescently-labelled anti-human IgG or IgM antibody and this three-part complex is visualised with the aid of a fluorescence microscope. To prevent interference from rheumatoid factors (Rf) and to reduce IgG competition in the IgM test, samples should be pre-treated with a suitable adsorbent reagent.

Kit Components

Materials Provided

1. Slides – 6 x10 well slides

SLIDE

On which recombinant Parvovirus B19 antigen expressed in insect cells has been stabilised. Foil pouch.

2. IgM Positive Control**

CONTROL + IgM

1 x 250µl positive sera containing specific anti-Parvovirus B19 antibodies. Contains thiomersal (0.01%). Green Cap.

3. IgG Positive Control**

CONTROL + IgG

1 x 250µl positive sera containing specific anti-Parvovirus B19 antibodies. Contains thiomersal (0.01%). Blue cap

4. Negative Control**

CONTROL -

1 x 250µl negative control sera containing no detectable anti-Parvovirus B19 antibodies. Contains thiomersal (0.01%).

Red cap

5. Fluorescent anti-IgM Conjugate

CONJ IgM

1 x 2ml Rabbit anti-human IgM antibody conjugated to fluorescent isothiocyanate (FITC). Contains Evans Blue counterstain and thiomersal (0.01%).

Purple cap

6. Fluorescent anti-IgG Conjugate

CONJ IgG

1 x 2ml Rabbit anti-human IgG antibody conjugated to fluorescent isothiocyanate (FITC). Contains Evans Blue counterstain and thiomersal (0.01%).

Yellow cap.

7. Wash Buffer Concentrate

BUF WASH 25X

1 x 45ml of concentrated (25x) Tris buffered saline with Tween 20 (0.25%) and thiomersal (0.01%).

Clear cap

8. Mounting Media

MM

1 x 2ml Tris – glycerol buffer. Contains thiomersal (0.001%).

Orange cap

9. Product Insert

INS

Instructions for use.

** Potentially Biohazardous Material

Additional Materials Required

- High quality distilled or deionised water.
- Accurate pipettes, micropipettes and

disposable tips to deliver 5µl to 50µl, 50µl to 200µl.

- Serum collection equipment
- Timer
- Clean volumetric laboratory glassware.
- A suitable Rf and IgG adsorbent reagent with a minimum IgG binding capacity of 18mg/ml.
- 1L beaker
- Dilution tubes and minifuge tubes (0.5ml).
- Benchtop minifuge.
- Incubation tray containing moistened tissue paper.
- Wash bottles and wash tray.
- 37°C incubator.
- No. 1 coverslips (22 x 50mm).
- Fluorescence microscope with appropriate filter combination for FITC (excitation filter 495nm, barrier filter 515nm), a halogen light source is recommended.
- Wax pencil
- Fan

Precautions

Safety

- For *in vitro* use only.
- The kit is intended for use by qualified laboratory staff only.
- All reagents derived from human origin are considered POTENTIALLY BIOHAZARDOUS MATERIAL. The Positive and Negative Control sera were tested by FDA-cleared methods for HBsAg and antibodies to HIV 1/2 and HCV and found to be negative. However, since no test can provide complete assurance of the absence of virus, treat all controls as potentially infectious.
- Some reagents contain thiomersal which may be toxic if ingested.
- Avoid contact with Evans Blue (in FITC IgM and IgG conjugates) as it is a potential carcinogen. If skin contact occurs, flush with large volumes of water.

- Dispose of all clinical specimens, infected or potentially infected material in accordance with good laboratory practice. All such materials should be handled and disposed of as though potentially infectious.
- Residues of chemicals, preparations and kit components are generally considered as hazardous waste. All such materials should be handled and disposed of as though potentially infectious.
- Wear protective clothing, disposable latex gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- Do not pipette materials by the mouth and never eat or drink at the laboratory workbench.

Procedural

- Do not use kit or individual reagents past their expiry date.

- Do not mix or substitute reagents from different kit lot numbers.
- Deviation from the protocol provided may cause erroneous results.
- Allow all reagents to come to room temperature (20 – 25°C) and mix well prior to use.
- Avoid leaving reagents in direct sunlight and/or outside 2-8°C for extended periods.
- When staining multiple samples on a slide avoid cross contamination between samples by marking between wells with a wax pencil.
- Application of excess mounting media may cause blurred fluorescence.
- Always use clean, preferably disposable, glassware for all reagent preparation.
- Care must be taken not to contaminate components and always use fresh pipette tips for each sample and component.

- Do not scratch the well with the pipette tip or dropper.
- Before commencing the assay an identification and distribution plan should be established.

Storage and Stability

- The kit is stable until the expiry date indicated on the outer box label provided it is stored between 2 – 8°C.
- Diluted Wash Buffer is stable for 6 months when stored at 2 – 8°C. Discard immediately if microbial contamination occurs.
- All unused components should be returned to 2-8°C storage immediately after use.

Specimen Collection and Storage

Serum samples should be obtained using aseptic laboratory techniques. These samples can be stored for up to 24 hours at 2-8°C and at –20°C or lower for longer periods. Fresh or fresh-frozen serum samples may be used, however, care should be taken to avoid repeated freeze-thawing.

Specimen and Reagent Preparation

Reagent Preparation

The Wash Buffer is supplied as liquid concentrate (25x). Prepare Wash Buffer for use by diluting 40ml of Wash Buffer Concentrate in 960ml of distilled water. Store in a clean closed container at 2-8°C for up to 6 months.

All remaining reagents are supplied ready to use and are at working dilution.

Specimen Preparation and Adsorption

- Qualitative Test:

For an IgG test dilute the patient serum sample 1:64 in Wash Buffer.

For IgM tests, samples should be adsorbed with suitable adsorbent for IgG and Rheumatoid factor (with a minimum binding capacity of 18 mg/ml). A final serum dilution of 1:16 should be prepared before use.

- Semi Quantitative Test:

For IgG tests, the sample 'titre' can be determined by serially diluting the patient sample 1:64, 1:256, 1:1024 in wash buffer or until a '+ 1' grade of fluorescence is achieved (see "Interpretation of Results").

For IgM tests, the sample 'titre' can be determined by serially diluting the 1:16 adsorbed patient sample to 1:64, 1:256 in Wash

Buffer or until a '+1' grade of fluorescence is achieved (see "Interpretation of Results").

Note: The 'titre' of the sample is the dilution used to achieve a '+1' grade of fluorescence

Note: Diluted samples should not be stored. If a repeat test is needed a fresh preparation should be used.

Assay Procedure

1. Allow all components to equilibrate to room temperature (20-25°C) before use.
2. Remove the desired number of slides from their sachets.
3. Mark between wells with a wax pencil to avoid cross contamination.

4. Dispense 20µl of each diluted sample and 20µl of the ready to use positive and negative controls onto numbered wells.

5. Incubate slides in an incubation tray for 3 hours at 35-39°C and 100% humidity.

6. Rinse slides briefly along the edge with Wash Buffer using a wash bottle. Place slides in a wash tray containing Wash Buffer (20ml/slide) for 10 minutes at room temperature (20-25°C). Rinse with Wash Buffer from wash bottle and fan dry slides.

7. For an IgM test dispense one drop (40µl) of anti-IgM FITC conjugate to each well. For an IgG test dispense one drop (40µl) of anti-IgG FITC conjugate to each well.

8. Incubate slides in the incubation tray in the dark for 30 minutes at 35-39°C and 100% humidity.

9. Repeat the washing procedure described in step 6 above.

10. Add one small drop of mounting media to the centre of each well and apply a coverslip.

11. Examine under a fluorescence microscope using 400 x magnification.

Interpretation of results

Negative: A sample is considered to be negative for Parvovirus B19 IgM and IgG antibodies if there are no visible fluorescence of the VP1 aggregates ('Bunch of Grapes' morphology).

Positive: A sample can be considered positive for Parvovirus B19 IgM and IgG antibodies if a positive fluorescent result is obtained at a dilution of $\geq 1:16$ and $\geq 1:64$, respectively. Positive fluorescence is indicated by staining of the distinct VP1 protein aggregates ('Bunch of Grapes' morphology) and can be graded as follows:

+4 = Brilliant green fluorescence indicating very high titre Parvovirus B19 antibody response.

+3 = Bright green fluorescence indicating high titre Parvovirus B19 antibody response.

+2 = Green fluorescence indicating medium titre Parvovirus B19 antibody response.

+1 = Dull green fluorescence indicating weak titre Parvovirus B19 antibody response. This also indicates the end point dilution or 'titre' of the sample.

+/- = Faint green fluorescence indicating possible Parvovirus B19 antibody response – confirm using available clinical information and other detection methods.

Cells which do not contain VP1 protein are included in each well to allow comparison of positive and negative cells. These stain blood red with Evans Blue counterstain.

Aspecific: Some individuals may exhibit an aspecific IgM reaction which is shown by a weak

fluorescence with a poorly defined staining pattern at a 1:16 serum dilution. This result should be verified by an alternative method (e.g. EIA, Immunoblot).

Quality Control Criteria

The Calibrator (Positive Control) and Negative Control must always be included to determine the validity of test results. Results of an assay are considered valid if the following criteria are met.

1. The IgG and/or IgM Positive Control yields a fluorescence greater than or equal to +2.
2. The Negative Control yields no visible fluorescence of the VP1 aggregates ('Bunch of Grapes' morphology)

If the above criteria are not met the assay is considered invalid and must be repeated.

Limitations of Use

- Results must be correlated with the patient's clinical and epidemiological profile and other clinical laboratory results in making the diagnosis of Parvovirus B19 infection.
- A negative result does not exclude the possibility of Parvovirus B19 infection. The development of a detectable antibody response may occur some days after infection. In the case of suspected Parvovirus B19 infection a negative result should be followed up with a repeat test two weeks later.
- Test performance may be affected by deviation from the procedure, interpretation or recommended precautions.
- Test results of specimens from immunocompromised patients may be difficult to interpret.

Performance Characteristics

Clinical Evaluation

The following results were obtained when sera from patients with symptoms of recent B19 infection (seroconversion panels) were tested using the Biotrin IgM IFA and a RIA method.¹⁷

Table 1 Clinical evaluation results

Days after onset	Biotrin IgM IFA		IgM RIA	
	No. of IgM Positives	%IgM Positive	No. of IgM Positives	%IgM Positive
1 - 5	26/29	89.6	26/29	89.6
6 - 10	22/22	100	22/22	100
11 - 20	8/8	100	8/8	100
21 - 30	8/8	100	8/8	100
31 - 50	4/4	100	4/4	100
>50	2/9	22.2	NT	NT

NT = Not Tested.

Seroprevalence

A total of 201 samples from randomly selected healthy blood donors were tested using the Biotrin Parvovirus B19 Immunofluorescent Assay for Parvovirus B19 IgG and IgM antibodies (IgG=201, IgM=200). 138 serum samples were IgG positive, representing 68.7% of samples screened, while only 2 samples (1%) were IgM positive. These findings are in agreement with seroprevalence studies reported in the literature.^{20, 21}

Sensitivity and Specificity

STUDY 1

The Centres for Disease Control and Prevention (CDC), maintain a panel of Parvovirus B19 sera which serves as an internal control and reference for their testing program for Parvovirus B19. The samples were characterised by clinical symptoms associated with a Parvovirus B19 infection. The majority of these patients were women with

household exposures to symptomatic children, or teachers at elementary schools with outbreaks of erythema infectiosum. This panel has also been tested using the CDC Parvovirus B19 native virus EIA²², Biotrin's Parvovirus B19 IgG EIA and Biotrin's Parvovirus B19 IgM EIA. A selection of sera from the well characterised CDC panel (n=30) were tested using the Biotrin Parvovirus B19 IgG and IgM IFA. A comparison of the results from the different methods is provided in Tables 2 & 3.

Table 2 Comparison of Biotrin IgG & IgM IFA with CDC native virus assay.

		CDC native virus assay*							
		Parvovirus B19 IgG				Parvovirus B19 IgM			
		Pos	Neg	Equiv	Tot	Pos	Neg	Equiv	Tot
Biotrin Parvovirus B19 IgG & IgM IFA	Pos	22	0	0	22	12	0	0	12
	Neg	0	5	0	5	0	15	0	15
	Equiv	0	0	0	0	0	0	0	0
	Tot	22	5	0	27	12	15	0	27

Only 27 of the 30 samples were tested on the CDC native virus assay.

Table 3 Comparison of Biotrin IgG & IgM IFA with Biotrin Parvovirus B19 IgG & IgM EIA assays.

		Biotrin's Parvovirus B19 EIA							
		Parvovirus B19 IgG				Parvovirus B19 IgM			
		Pos	Neg	Equiv	Tot	Pos	Neg	Equiv	Tot
Biotrin Parvovirus B19 IgG & IgM IFA	Pos	22	0	0	22	12	0	0	12
	Neg	0	8	0	8	0	18	0	18
	Equiv	0	0	0	0	0	0	0	0
	Tot	22	8	0	30	12	18	0	30

Pos = positive

Neg = negative

Equiv = equivocal

Sensitivity = True Positives divided by (True Positives + False Negatives + Equivocals) x 100

% Sensitivity = 100%

Specificity = True Negatives divided by (True Negatives + False Positives + Equivocals) x 100

% Specificity = 100%

STUDY 2

76 samples were evaluated for the presence of IgG & IgM antibodies using 4 separate screening methods – Hillcrest Biologicals Parvovirus IgG & IgM EIA assay (EIA), Microbiology Reference Lab Western Blot assay (WB), VRL Colindale RIA assay (RIA) and Biotrin IFA assay (IFA). The 4 screening methods for each antibody (IgG & IgM) were compared. An average result was considered when the screening of at least 3 out of 4 methods agreed or when at least 2 out of 3 methods agreed (when the other method was equivocal).

For Parvovirus B19 IgG antibodies, 53 samples were determined to be positive and 23 samples negative. The Biotrin IgG IFA correctly identified the 53 positive samples, yielding a sensitivity of 100%, however the IFA incorrectly identified 2 of the negative samples yielding a specificity of 92%.

For Parvovirus B19 IgM antibodies, 26 samples were determined to be positive and 50 samples negative. The Biotrin IgM IFA correctly identified the 50 negative samples, yielding a specificity of 100%, however the IFA incorrectly identified 1 of the positive samples yielding a sensitivity of 96%.

Cross Reactivity

100 serum samples were screened in order to establish the specificity of the Biotrin IgM IFA. All the samples were obtained from patients with the following diseases:

Table 4 Cross-reactive studies

Virus	Number of Positives
Varicella Zoster Virus	0/5
Rheumatoid Factor	0/5
Mumps	0/5
Rubella	0/9
Rubeola	0/5
CMV	0/13
EBV	0/12
Herpes Simplex virus - 1	0/5
Herpes Simplex virus - 2	0/5
Toxoplasma gondii	0/11
Lyme	0/5
Thyroiditis	0/5
Lupus	0/5
Rheumatoid Arthritis (RA)	0/5
ANA	1/5*

* This sample also produced a positive result using Biotrin's Parvovirus B19 IgM EIA & Parvoblot IgM Westernblot assay.

Reproducibility

• Intra – assay reproducibility

The following results were obtained when negative, low and high titre IgG and IgM samples were screened in three kits from one production lot.

Table 5 Intra-assay reproducibility

DEGREE OF FLUORESCENCE

	Slide No.	IgG Positive	IgG Cut Off	Negative	IgM Cut Off	IgM Positive
Kit 1	1	+4	+1/2	-	+1	+2/3
	2	+4	+1/2	-	+1	+2/3
	3	+4	+1/2	-	+1	+2/3
	4	+4	+1/2	-	+1	+2/3
	5	+4	+1/2	-	+1	+2/3
	6	+4	+1/2	-	+1	+2/3
Kit 2	1	+4	+1/2	-	+1	+2
	2	+4	+2	-	+/-	+2
	3	+4	+1/2	-	+/-	+2
	4	+4	+1/2	-	+1	+2
	5	+4	+1/2	-	+/-	+2
	6	+4	+1/2	-	+/-	+2
Kit 3	1	+3/4	+1/2	-	+/-	+2
	2	+3	+1	-	+/-	+2
	3	+3	+1	-	+1	+2
	4	+3/4	+1/2	-	+/-	+2
	5	+3	+1/2	-	+/-	+2
	6	+3	+1	-	+1	+2

• **Inter – site reproducibility**

The following results were obtained when negative, low, medium and high titre IgG (Table 6) and IgM (Table 7) samples were screened in triplicate on 3 separate days at two independent evaluation sites

Table 6 Inter – site reproducibility (IgG)

IgG Sample	Site 1			Site 2		
	3.11.94	8.11.94	9.11.94	10.10.94	11.10.94	12.10.94
1	neg	neg	neg	neg	neg	neg
1	neg	neg	neg	neg	neg	neg
1	neg	neg	neg	neg	neg	neg
2	neg	neg	neg	neg	neg	neg
2	neg	neg	neg	neg	neg	neg
2	neg	neg	neg	neg	neg	neg
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
4	pos	pos	pos	pos	pos	pos
4	pos	pos	pos	pos	pos	pos
4	pos	pos	pos	pos	pos	pos
5	Wk pos	Wk pos	Wk pos	pos	pos	pos
5	Wk pos	Wk pos	Wk pos	pos	Wk pos	pos
5	Wk pos	Wk pos	Wk pos	pos	pos	pos
6	pos	Wk pos	Wk pos	Wk pos	Wk pos	Wk pos
6	pos	Wk pos	Wk pos	Wk pos	Wk pos	Wk pos
6	pos	Wk pos	Wk pos	Wk pos	Wk pos	Wk pos
7	Wk pos	pos	pos	pos	pos	pos
7	pos	pos	pos	pos	pos	pos
7	pos	pos	pos	pos	pos	pos
Positive Control	pos	pos	pos	pos	pos	pos
Negative Control	neg	neg	neg	neg	neg	neg

Wk Pos =weak positive

Pos = positive

Neg = negative

Table 7 Inter – site reproducibility (IgM)

IgM Sample	Site 1			Site 2		
	3.11.94	8.11.94	9.11.94	10.10.94	11.10.94	12.10.94
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	Wk pos
2	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
4	Wk pos	Wk pos	pos	Wk pos	pos	pos
4	Wk pos	Wk pos	pos	Wk pos	pos	Wk pos
4	Wk pos	Wk pos	pos	Wk pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
7	neg	Wk pos	Wk pos	neg	neg	neg
7	neg	Wk pos	Wk pos	neg	neg	neg
7	neg	Wk pos	Wk pos	neg	Wk pos	neg
Positive Control	pos	pos	pos	pos	pos	pos
Negative Control	neg	neg	neg	neg	neg	neg

Wk Pos = weak positive

Pos = positive

Neg = negative

Summary of Parvovirus B19 IF IgG & IgM Procedure

Important Note:

Please read the entire product instruction leaflet before starting the assay.

This summary is for quick reference only.

- For an IgG test dilute patient serum samples 1 in 64 in Wash Buffer
- For an IgM test adsorbed samples should be diluted to final dilution of 1 in 16
- Pipette 20µl of Positive Control (ready to use), 20µl Negative Control (ready to use) and 20µl of prepared samples onto wells
- Incubate for 3 hours @ 35-39°C, 100% humidity
- Wash slides with Wash Buffer
- For IgG test add 40µl IgG conjugate to each well

- For IgM test add 40µl IgM conjugate to each well
- Incubate for 30 min in the dark @ 35-39°C, 100% humidity
- Wash slides with Wash Buffer
- Add mounting media to each well
- Examine under a fluorescence microscope

Interpretation of Symbols

	- In vitro diagnostic medical device		- Use by end of
	- Temperature limitation		- Batch code
	- Manufacturer		- Catalogue Number

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Campione IgM	Sede 1			Sede 2		
	3.11.94	8.11.94	9.11.94	10.10.94	11.10.94	12.10.94
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	Deb pos
2	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
4	Deb pos	Deb pos	pos	Deb pos	pos	pos
4	Deb pos	Deb pos	pos	Deb pos	pos	Deb pos
4	Deb pos	Deb pos	pos	Deb pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
7	neg	Deb pos	Deb pos	neg	neg	neg
7	neg	Deb pos	Deb pos	neg	neg	neg
7	neg	Deb pos	Deb pos	neg	Deb pos	neg
Controllo positivo	pos	pos	pos	pos	pos	pos
Controllo negativo	neg	neg	neg	neg	neg	neg

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Amostra IgM	Centro 1			Centro 2		
	3.11.94	8.11.94	9.11.94	10.10.94	11.10.94	12.10.94
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
1	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	Frac Pos
2	pos	pos	pos	pos	pos	pos
2	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
3	pos	pos	pos	pos	pos	pos
4	Frac Pos	Frac Pos	pos	Frac Pos	pos	pos
4	Frac Pos	Frac Pos	pos	Frac Pos	pos	Frac Pos
4	Frac Pos	Frac Pos	pos	Frac Pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
5	pos	pos	pos	pos	pos	pos
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
6	neg	neg	neg	neg	neg	neg
7	neg	Frac Pos	Frac Pos	neg	neg	neg
7	neg	Frac Pos	Frac Pos	neg	neg	neg
7	neg	Frac Pos	Frac Pos	neg	Frac Pos	neg
Controlo positivo	pos	pos	pos	pos	pos	pos
Controlo negativo	neg	neg	neg	neg	neg	neg

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9. Jordan JA.
Identification of human Parvovirus B19 in idiopathic nonimmune hydrops fetalis.
Am J Obstet Gynecol 174: 37 – 42, 1996.
10. Centres for Disease Control.
Risks Associated with Human Parvovirus B19 Infection.
MMWR CDC Surveill. Summ., 38: 81 – 89, 1989.
11. Woernle CH, Anderson LJ, Tattersall P, et al.:
Human Parvovirus B19 during pregnancy.
J. Infect. Dis. 156: 17-20, 1987.
12. Gratacos E, Torres PJ, Vidal J, et al.:
The incidence of Human Parvovirus B19 infection during pregnancy and its impact on perinatal outcome.
J Infect. Dis. 171: 1360 – 1363, 1995.
13. Public Health Laboratory Service Working Party on Fifth Disease.
Prospective Study of Human Parvovirus (B19) infection in pregnancy.
Br. Med. J. 300: 1166 – 1170, 1990.
14. Rodis JF, Hovick TJ, Quinn DL, et al.:
Human Parvovirus Infection in pregnancy.
Obstet. Gynecol. 72: 733-738, 1988.
15. Brown KE, Hibbs JR, Gallinella G, et al.:
Resistance to Parvovirus B19 infection due to lack of virus receptor (Erythrocyte P Antigen).
N. Eng. J. Med. 330: 1192-1196, 1994.
16. Anderson MJ, Khoussam MN, Maxwell DJ, et al.:
Human Parvovirus B19 and Hydrops Fetalis.
Lancet 1: 535, 1988.
17. Coons AH, Kaplan MM, et al.:
J. Exp. Med. 91: 1-13, 1950.
18. Brown CS, van Bussel MJ, Wassenaar AL, et al.:
An immunofluorescence assay for the detection of Parvovirus B19 IgG and IgM antibodies based on recombinant viral antigen.
J. of Virological Methods 29: 53-62, 1990.
19. Browns CS, Salimans MM, Noteborn MH, et al.:
Antigenic Parvovirus B19 coat proteins VP1 and VP2 produced in large quantities in a Baculovirus expression system.
Virus Res. 51: 197-212, 1990.
20. Searle K, Gulliard C, Enders G, et al.:
Parvovirus B19 Diagnosis in Pregnant Women - Quantification of IgG Antibody Levels (IU/ml) with Reference to the International Parvovirus B19 Standard Serum.
Infection 25: 32-34, 1997.
21. Gay NJ, Hesketh LM, Cohen BJ, et al.:
Age specific antibody prevalence to Parvovirus B19: how many women are infected in pregnancy?
Comum. Dis. Rep. CDR Rev. 4(9): R104-107, 1994.
22. Anderson LJ, Tso C, Parker RA, et al.:
Detection of antibodies and antigens of human Parvovirus B19 by enzyme-linked immunosorbent assay.
J. Clin. Microbiol. 24(4): 522-526, 1986.