



Please note the following information relevant to the detection of the H1N1 (Swine Flu) antigen by the Biotrin [Respiratory Virus Panel IFA kit](#).

**In the [Biotrin RVP kit \(code V4RVP\)](#) it has been verified that the monoclonal antibodies directed against influenza A have positively detected the new strain of H1N1, most recently known as "swine flu". The new strain was shown to have been detected in both direct specimens and after amplification in culture.**

**As per the limitations of use in the V4RVP product information for use booklet it is advised that a negative result does not preclude the presence of influenza A. Under the same limitation therefore a negative result does not exclude infection with swine-origin influenza A (H1N1) virus. Results should be interpreted in conjunction with patient related clinical information and other diagnostic procedures.**

***Important Note for United States Sales and Usage Limitation: For Research Use Only. Not for use in diagnostic procedures.***

Finally please do not hesitate to [contact us](#) if you have any questions,

Thank you